EXCEL ROADSTER DX8 DELUXE Specification Sheet



Full Specifications

- Seat Width46 cmSeat Depth40 cmSeat Height53 58 cm*(From Ground)132 cmTotal Length132 cmTotal Width61 cmTotal Height102 cm
 - Maximum Speed Range Turning Radius Ground Clearance Battery Charger Total Weight (Including Battery) Maximum User Weight

8 mph 43 km (27 miles)** 153 cm 11 cm 2 x 12V / 55Ah 24V / 4Ah 92 kg

136 kg (21 Stone)

KEY: * = Adjustable or optional ** = Varies depending on weather and terrain conditions

www.vanosmedical.co.uk