

DRIFFIELD MOBILITY



Welcome To Your New Pavement Scooter



Here are some handy tips to help you get the most out of your pavement scooter.



DO

- Remember your pavement scooter doesn't like getting wet.
- Keep to pavements, your Class 2 scooter cannot be driven along the road.
- Store your scooter somewhere dry, either indoors, garage or shed.
- Turn the power off when transferring to and from the scooter.
- Avoid rough or soft terrain and long grass wherever possible.
- Turn the power off when not in use.
- Keep your battery regularly charged for battery health.
- Buy a scooter cover to protect from the elements if regularly used outdoors.
- Remember to read the manufacturers instructions.



DON'T

- Leave your scooter outside in the rain or use during rainfall without a cover.
- Exit or enter the scooter unless the power is switched off.
- Switch off the power whilst driving – this will cause the scooter to stop abruptly.
- Exceed the stated weight limit.
- Travel through water or over soft terrain (mud, sand etc.)
- Turn sharply at full speed.
- Leave on charge for more than 12 hours .

Thank you for your purchase .