

DRIFFIELD MOBILITY



Welcome To Your New Road Scooter



Here are some handy tips to help you get the most out of your Road scooter.



DO

- Store your scooter somewhere dry, either indoors, garage or shed.
- Be seen at night. Always use your scooter lights for high visibility
- Turn the power off when transferring to and from the scooter.
- Avoid rough or soft terrain and long grass wherever possible.
- Turn the power off when not in use.
- Keep your battery regularly charged for battery health.
- Buy a scooter cover to protect from the elements.
- Register your new road scooter with DVLA
- Remember to read the manufacturers instructions.
- Get your road scooter serviced annually.



DON'T

- Ride along the pavement unless at 4mph or below.
- Exit or enter the scooter unless the power is switched off.
- Switch off the power whilst driving – this will cause the scooter to stop abruptly.
- Exceed the stated weight limit.
- Travel through water or over soft terrain (mud, sand etc.)
- Turn sharply at full speed.
- Leave on charge for more than 12 hours .

Thank you for your purchase .