

New!

Get Out There!

GO-CHAIR

GET OUT THERE IN STYLE WITH THE ALL-NEW GO-CHAIR



Shown in: Ruby Red

The all-new Go-Chair is re-engineered from the ground up, offering a bold new style available in an array of contemporary colors. Enhanced performance and comfort, along with feather-touch disassembly, allows you to enjoy light-weight travel and independence on the go. With an increased weight capacity of 300 lbs. and a sleek, bold look, the new Go-Chair makes travel easy.



One-hand "Feather Touch" disassembly



Variety of color options

SPECIFICATIONS

GO CHAIR


Drive Wheels	9" (solid)
Front Casters	6" (solid)
Rear Anti-tips	3" (solid)
Maximum Speed ¹	Up to 3.7 mph
Ground Clearance ²	1.6" at anti-tip bracket 1.9" at frame
Turning Radius ²	25.4" with foot platform
Length	33" with front rigging
Width ⁶	22.6"
Seat Type	Foldable, compact, lightweight
Seat to Ground Height	22"-23"
Batteries ^{3,5}	(2) 18 Ah, 12 V, deep-cycle
Range ^{1,4}	Up to 8.74 miles
Charger	Off-board, 2 amp
Electronics	40 amp, PG GC3 (non-programmable)
Weight Capacity	300 lbs.
Base Weight	67 lbs. (front & rear)
Standard Seat Weight	30 lbs.
Battery Pack Weight ⁵	31 lbs.
Suspension	Limited
Heaviest Piece	36 lbs. (rear)

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition tire type and condition. These specifications can be subject to a variance of +10%, -5%.
² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
³ AGM or Gel-Cell type required. *See "Batteries and Charging."
⁴ Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
⁵ Battery weight may vary based on manufacturer.
⁶ Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA and ISO 7176 series standards. This product is not coded by DME PDAC or does not meet criteria: K0899. All specifications subject to change without notice.

COMPATIBLE LIFTS

- Outlander Full Platform
- Olympian
- Backpacker Series
- Commander

 To find a compatible lift for your vehicle/mobility product, use our new tool: pridemobility.com/lift-finder/

4-wheel

FEATURES

- Choice of 6 colors
- 18" x 17", 5-point swivel seat with latch release
- Increased 300 lb. weight capacity
- Faster speed - up to 3.7 mph
- Larger, 18 amp batteries
- Improved agility - longer wheel base; wider stance; compact turning radius 25.4"
- Dual swing-away, under seat storage bins
- Larger, height & width adjustable armrests
- Larger footplate
- One hand feather-touch disassembly
- 60" lap belt

OPTIONS

- Cup Holder
- Lap Belt- 50", 70"
- Safety Flag
- Oxygen Holder
- Dual Saddle Bag
- Weather Cover, Micro
- Additional Battery Pack
- Rear Mount Bracket
Cane/Crutch Holder
Single (3"),
Double (3")
- Armrest Mount Cane Clip

Pride
Mobility Products Co.
— LIVE YOUR BEST! —

GO GO
TRAVEL MOBILITY

Beamsville, Ontario • 888-570-1113

[PrideMobility](#) [PrideMobilityProductsCorp](#) [PrideMobilityProductsCanada](#)

pridemobility.com