

**ROMA**  
MEDICAL



P319 - Reno II Powerchair

**ROMA**  
MEDICAL

# ROMA

# RENO II Powerchair



Detachable battery pack



PG Drives controller



Storage basket and flip up footrest



Easily dismantled without tools

## Special Features

- UK designed & registered by Roma.
- New battery pack 2 x 22ah.
- New patented docking system.
- Intelligent, regenerative, electromagnetic brakes.
- New frame locking mechanism.

Available in **Blue**



## Specifications

<b>Model No.* -</b>	<b>P319</b>	<b>Safe Working Slope -</b>	6 degrees / 10%
Class (Invalid Carriage)	2	<b>Maximum User Weight -</b>	136kgs / 21 stone
Overall Length -	869mm / 35.2"	<b>Battery Size 12v x 2 -</b>	22Ah
Overall Width -	571mm / 22.5"	<b>Range** -</b>	20km / 12.5 miles
Overall Height -	860mm / 34"	<b>Motor -</b>	DC24V, 90W
Overall Height (seat removed) -	427mm / 16.8"	<b>Maximum Speed -</b>	4mph / 6.4 kmph
Seat Width -	460mm / 18"	<b>Wheel / Tyre Size (Front) -</b>	150mm / 6"
Seat Depth -	380mm / 15"	<b>(Back) -</b>	230mm / 9"
Seat Height (to ground) -	532mm / 21"	<b>Turning Radius -</b>	701mm / 27.6"
Seat Height (to floorpan) -	427mm / 16.8"	<b>Battery Charger -</b>	4A off board
Ground Clearance -	52mm / 2"	<b>Battery Pack Weight -</b>	14kgs / 31 lbs

\* Lower limb amputees should refer to user manual before use. \*\* Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel. \*\*\* Measurement to underside of chassis



# ROMA MEDICAL

York Road Bridgend Industrial Estate Bridgend CF31 3TB  
Telephone : 01656 674488 Facsimile : 01656 674499  
Email : sales@romamedical.co.uk

## www.romamedical.co.uk

Roma Medical reserve the right to change specifications without prior notice

Dealer Stamp

