



# Denver & Denver **PLUS**

Model S130 & S135 | Class 2 Scooter

Available in:   

# Denver & Denver <sup>PLUS</sup>

Model S130 & S135 | Class 2 Scooter



Available in Red, Orange or Graphite Grey.



Front and rear shock absorbers.



Front and rear mud guards.



Shopping basket.



Detachable battery pack.



Easy to use controls.



York Road  
Bridgend Industrial Estate  
Bridgend  
CF31 3TB

Tel: 01656 674488  
Fax: 01656 674499  
Email: [sales@romamedical.co.uk](mailto:sales@romamedical.co.uk)

[www.romamedical.co.uk](http://www.romamedical.co.uk)

 We are social. Follow us on Facebook  
[/RomaMedicalAids](https://www.facebook.com/RomaMedicalAids)



\*Lower limb amputees should refer to the user manual before use.  
\*\* Passenger weight, ground conditions and battery condition/size can affect distance of travel.  
\*\*\* Measurement of underside of chassis.

Roma Medical reserve the right to change specifications without prior notice.

## Special Features

Front LED headlight.

Front and rear shock absorbers.

Front and rear mud guards.

Quickly dismantles without tools for easy transport and storage.

Anti-tipping wheels at rear.

Comfortable swivel adjustable height seat with rear storage pocket.

Flexible seat slider.

Flip up and width adjustable armrests.

Adjustable angle tiller.

Detachable battery pack.

Durable rubber foot mat and foot rests.

Two charging points, one on-board and one off-board.

Generous leg room.

Available in Bright Orange, Bright Red and Graphite Grey.

## Specifications (General)

Class (Invalid Carriage):	2
Overall Length:	1120mm / 44.1"
Overall Width:	550 - 650mm / 21.7 - 25.6"
Overall Height:	930 - 970mm / 36.7 - 38"
Seat Width:	440mm / 17"
Safe Working Slope:	6° / 10%
Maximum User Weight:	136 kg / 21.5 st
Maximum Speed:	6 kmph / 4 mph
Wheel Size Front:	7.5"
Wheel Size Rear:	8.5"
Turning Radius:	1440mm / 56.7"
Battery Charger:	2 Amp Standard
Ground Clearance***:	100mm
Colours:	Orange, Red & Graphite Grey

## Specifications (Denver)

Model No.*:	S130
Range*	16 km / 10 miles
Battery Size:	2x 12ah
Overall Weight	56 kg / 124 lbs

## Specifications (Denver Plus)

Model No.*:	S135
Range*	32 km / 20 miles
Battery Size:	2x 38ah
Overall Weight:	73 kg / 160 lbs