



NEW PRODUCT
UPDATE..

SLEEK STYLE AND
LUXURIOUS FEATURES



PRIDE® **MOBILITY** SCOOTERS
LIVE YOUR BEST®

Pride®
Mobility Products Ltd.
— LIVE YOUR BEST® —



Tel: 01869 324600

Email: sales@pride-mobility.co.uk

Web: www.pride-mobility.co.uk

(Version 2: June 2017)


Features and Specification

Introduction

The Apex Rapid with CTS Suspension (Comfort-Trac Suspension) delivers luxury and convenient features like a wraparound delta tiller, alloy style 8" front and rear tyres, high level charging port, LED kerb light, storage basket and a compact easy to transport design. The Apex Rapid's sophisticated design and abundant features make it an exceptional travel mobility choice.



CTS Front and Rear Suspension



Also Available
with a:
12 Ah Battery



Disassembles
into 5 pieces

4x

Other optional shroud
colours available.



Yellow



Black



Orange



Pink

Available Colours

The Apex Rapid comes in two optional colours:














Red



Blue

Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
										
21 Stone 6lbs. (137 kg)	4 Mph (6.4 km/h)	3.81 cm (1.5 inch's)	113.6 cm (44.7 inch's)	101 cm (39.7 inch's)	52.7 cm (20.5 inch's)	Front: 20cm Rear: 20cm	6.23 Stone 39.6k.g.	17.8 k.g. Front Section	Folded Seat 10.4k.g	x2 18 Ah Battery

Warranty

Two-year limited warranty on frame
One-and-a-half-year limited warranty on Controllers
One-year limited warranty Electronic Components,
battery and Charger.

Notes

Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
Battery weight may vary based on manufacturer.

Available on



Tel: 01869 324600

Email: sales@pride-mobility.co.uk

Web: www.pride-mobility.co.uk

Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST —