U-DRIVE

WHEELCHAIR TO POWERCHAIR CONVERSION MODULE





 ϵ



The unique U-Drive system converts the XS aluminium self propelled wheelchair, into a user controlled electric powerchair.

With U-Drive, U can propel the wheelchair or if required simply attach the module and U drive.

- The U-Drive can easily be fitted or removed.
- The excellent P G Drives S-Drive Joystick provides true user control.
- The S-Drive controller has a fully programmable & comprehensive fault management system.
- New push and slide battery removal system.
- Very portable, making transportation simple.
- When detached the U-Drive and wheel chair will fit into most car boots.
- Maximum speed of 4mph*.
- Maximum range of 6 miles*.
- Simple to manoeuvre and reverse.
- The U-Drive has a powerful 150 watt, 24V motor, capable of coping with most terrains.
- Maintenance free airline friendly 24V 12ah valve regulated sealed lead acid battery pack.
- On or off board charging.
- 2A charger as standard for off-board charging.
- Heaviest part 22kg.
- Dual control option available.

*Subject to conditions of use





EASILY FITTED OR DETACHED

SIMPLE BATTERY REMOVAL

CODE	TOTAL WEIGHT	WEIGHT WITHOUT BATTERY	WEIGHT BATTERY PACK ONLY	MAXIMUM GRADIENT	MAXIMUM RANGE	WEIGHT CAPACITY
U-Drive	30.5kg	22kg	8.5kg	6°	6 Miles	115kg (18st)

The unique U-Drive system converts the XS aluminium wheelchair, into a user controlled electric powerchair. This revolutionary new product is an excellent addition to the Powerstroll range from Drive Medical.

For more information please telephone 01422 314488 or Email: enquiries@drivemedical.co.uk

