## POWERSTROLL Wheelchair Accessories



"Takes the strain out of pushing a wheelchair!"

CE

This excellent Powerstroll will turn most manual wheelchairs into carer controlled powered wheelchairs. It takes the strain out of pushing a heavy wheelchair and makes easy work of hills, ramps and uneven surfaces. Its unique one click docking station (patented) means that it can be fitted or removed easily in seconds.



The Powerstroll can easily be fitted and removed within seconds through its unique one-click docking station and the new push and slide battery removal system.



Easy to fit & remove





Hand control with variable speed dial, forward/reverse function, illuminated battery gauge and fingertip trigger control.



Fits most standard wheelchairs up to 20" seat width

## Superb Comfort, Style & Performance...

- » Allows speed of up to 4mph\*
- » Travels up to 10 miles (16 km) on a full battery charge\*
- » The Powerstroll can easily be fitted or removed within seconds
- » The excellent P G Drives S-Drive controller with fully programmability & comprehensive fault management system
- » New push and slide battery removal system
- » It is very portable, making transportation simple
- » Fits most standard wheelchairs up to 20" seat width
- » Additional parts available for wider wheelchairs
- Hand control with variable speed dial, forward/reverse function, illuminated battery gauge and fingertip trigger control

- Unique one-click docking station for ease of assembly (patented)
- Comes complete with a Carry Holdall with shoulder strap
- » Simple to manoeuvre and reverse
- The Powerstroll has a powerful 150 watt (200 watt-HD), 24V motor, capable of coping with most terrains
- Maintenance free airline friendly 24V 12ah valve regulated sealed lead acid battery pack
- » 2A charger as standard for off-board charging
- » Weight capacity 135kg (21st), this includes the wheelchair and user combined
- Heavy Duty model weight capacity 180kg (28st)
- \* Subject to conditions of use

## CODE: PWCPPoo9 & PWCPPoo9HD DESCRIPTION: Powerstroll S Drive



## Leaflet Code: LL144

If you would like to view this product please contact your local stockist, details below.