

**Z23031 - HEIGHT ADAPTOR BRACKET FOR PWCPP010**

These brackets will allow the telescopic bar at the front of the powerpack to be installed at its standard height or at a new height approx 2" (5cm) lower.

This is ideal for installation on wheelchairs have a high lower horizontal bar or lack a hemi-shape (such as the XS2 and SD2).

The brackets are supplied in pairs.
No further fasteners are required.

Z23032 - HEIGHT ADAPTOR BRACKET FOR PWCPP009 AND PWCPP009HD

These brackets will allow the telescopic bar at the front of the S-Drive powerpack to be installed at its a new height approx 2" (5cm) lower than the standard height,

This is ideal for installation on wheelchairs have a high lower horizontal bar or lack a hemi-shape (such as the XS2 and SD2).

The brackets are supplied in pairs with all required fasteners.

**Z23017 - MOUNTING BRACKETS FOR SENTRA**

These brackets allow the PWCPP009 and PWCPP009HD to be installed on the Sentra wheelchair. The brackets are at a fixed height tailored to the Sentra.

The brackets are supplied in pairs with all required fasteners.

Z23014-08 - EXTENSION BAR

Extension bar to extend the telescopic bar on the PWCPP010, PWCPP009 and PWCPP009HD so it will fit wheelchairs with up to a 24" seat width.

