



# Reno Elite

Model P325 & P325/CS | Class 2 Powerchair

Available in:  

# Reno Elite

Model P325 | Class 2 Powerchair



Available in Flame Red or Flame Blue.



Luxury padded captain seat.



Rear mud guards.



Storage basket.



Dynamic shark controller.



Available in a standard seat variant.



York Road  
Bridgend Industrial Estate  
Bridgend  
CF31 3TB

Tel: 01656 674488  
Fax: 01656 674499  
Email: sales@romamedical.co.uk

[www.romamedical.co.uk](http://www.romamedical.co.uk)

 We are social. Follow us on Facebook  
/RomaMedicalAids



\*Lower limb amputees should refer to the user manual before use.  
\*\* Passenger weight, ground conditions and battery condition/size can affect distance of travel.  
\*\*\* Measurement of underside of chassis.

Roma Medical reserve the right to change specifications without prior notice.

## Special Features

- Luxury padded Captain Seat.
- Easily dismantles without tools.
- Compact and lightweight design.
- Dynamic Shark controller.
- On and off-board charging.
- Intelligent, regenerative, electromagnetic brakes.
- Non-marking puncture proof tyres.
- Removable battery pack.
- Rear mud guards.
- Height adjustable, 180° swivel seat for easy transfer.
- Storage Basket.
- Flip-up armrests.
- Lapbelt as standard.
- Crash Tested.\*
- Also available in a standard seat variant.
- Available in Flame Red or Flame Blue.

\*Only the Captain Seat (P325/CS) variant has been crash tested.



## Specifications (General)

Overall Length:	880mm / 34.6"
Overall Width:	600mm / 23.6"
Max Height:	1000mm / 39.3"
Safe Working Slope:	6° / 10%
Maximum User Weight:	136kg / 21st
Battery Size:	22Ah
Range**:	20km / 21.5 miles
Maximum Speed:	6 kmph (4 mph)
Wheel Size (Front):	150mm / 6"
Wheel Size (Back):	250mm / 10"
Turning Radius:	550mm / 21.6"
Battery Charger:	2 Amp
Ground Clearance***:	50mm / 1.96"
Colours:	Flame Red & Flame Blue

## Specifications (Captain Seat)

Model No.*:	P325/CS
Seat Width:	460mm / 18.1"
Backrest Height:	490mm / 19.3"
Overall Weight:	67kg

## Specifications (Standard Seat)

Model No.*:	P325
Seat Width:	440mm / 17.3"
Backrest Height:	390mm / 15.3"
Overall Weight:	62kg